

# Privacy-Preserving Fall Detection System



Romain Labbe, Mathys Bitsch, Adrian Martinez Gao

# Falls: a Public Health Challenge

1/3

Adults aged 65+ fall each year


20k

Deaths from falls, France, 2024 (+18%)

Leading cause of accidental death among older adults

# A Growing Issue

Deaths from falls among adults aged 65 and over



Source: Santé Publique France <https://www.santepubliquefrance.fr/chute/donnees>

# | The "Long Lie" Fall – A Critical Emergency



**Delayed Rescue:** Lying on the floor for more than an hour drastically increases complications.



**Consequences :** Dehydration, muscle and tissue damage, hypothermia, delayed medical treatment, psychological trauma



**The 6-months toll :** **50% mortality rate within 6 months** following a "Long Lie" fall (>1h) for adults over 65, even without direct injuries from the fall

# The Privacy Concern



## **Cameras Opt-Out**

80% to 90% of older adults refuse optical cameras in their bedroom or bathroom.





## **High-Risk Area**

Bathroom is the #1 location for falls, yet it is also the most private space.



# Solution: mmWave Radars





# Technologies Comparison



| Metric               | Optical Cameras            | Wearables (Bracelets)      | Radar mmWave             |
|----------------------|----------------------------|----------------------------|--------------------------|
| Privacy Preservation | ✗ Non                      | ✓ High                     | ★ Superior (Anonymized)  |
| Auto-Detection       | ✓ Yes                      | ✗ Requires Manual Action   | ✓ Fully Automatic        |
| Low Light / Bathroom | ✗ Limited                  | ✓ Yes                      | ✓ Yes (RF Sensing)       |
| User Compliance      | ✗ Low (Privacy resistance) | ⚠ Medium (Often forgotten) | ✓ High (Passive sensing) |


# Multitracking

**Presentation**










# Fall Detection

**Presentation**



**Question Time**